

中科瑞泰(北京)生物科技有限公司

Tel: 010-58437227 82598075 Fax: 010-82597807 http:// <u>www.real-tims.com.cn</u> E-mail: <u>real-times@163.com</u>

Amphotericin B solubilized from Streptomyces sp

可溶性两性霉素 B

CAS Number: 1397-89-3 C₄₇H₇₃NO₁₇ MW 924.08

产品编号	粉末含量	两性霉素B有效含量
AE437-25mg	55mg	25mg

Brand: AmrescoE437

Storage: 2-8°C

• Specification:

Form.....fine, light yellow, lyophylized powder

Stock solution......2.5mg/ml(以有效含量计算) in H₂O

Product description:

Amphotericin B is a polyene antifungal antibiotic from *Streptomyces sp.* It has a high affinity for sterols, primarily ergosterols, of fungal and bacterial cell membranes. After binding to sterols, it forms channels in the membranes, causing small molecules to leak out. Amphotericin B is effective against fungi and yeast. Normal usage for maintenance of cell cultures is 2.5 mg/l with penicillin and streptomycin in the medium. For cultures already contaminated with yeast and fungus, it is recommended to use this product at 2-4 times the normal level (5-10 mg/l), without penicillin and streptomycin for 2-3 subcultures. Once the infection is under control, normal maintenance levels of amphotericin B should be used. Amphotericin B is not soluble in H₂O but in DMSO. Amphotericin B solubilized is a colloidal suspension of Amphotericin B, using sodium deoxycholate as the solubilizing agent. The product is approximately 45% amphotericin B, 35% sodium deoxycholate; the balance being sodium phosphate and sodium chloride.

• Preparation instructions:

Amphotericin B solubilized is reconstituted at the water in a 2.5 mg/ml stock solution. The solution must be filtered to sterilize by 0.22µm membrane.

• Storage/Stability:

The lyophylized powder amphotericin B solubilized can be stable for 5 years under -20 $^\circ\!C$ protected from air and light. Stock solutions at 2-8 $^\circ\!C$ are stable for up to 1 month.

Reference:

1. Bolard, J., Biochim. Biophys. Acta, 864: 257-304 (1986).

2. Venegas, B., et al., Biophys. J., 85, 2323-32 (2003).

- 3. Vertut-Doi, A., et al., Biochem. Biophys. Res. Commun., 157, 692-97 (1988).
- 4. Sau K., et al., J. Biol. Chem., 278, 37561-68 (2003).
- 5. Perlman, D., Methods Enzymol., 58, 110-16 (1979).